

Abstract

Detection and analysis of two or more co-localized antigens by light microscopy has been hindered by the difficulty in discerning and quantifying overlapping chromogens. Performing multi-analyte immunohistochemistry has many potential benefits and applications in the field of drug target evaluation and patient screening, including: 1) determination and quantification of immunoreactivity co-localization; 2) evaluation of spatial relationships between immunoreactive cell types; 3) signal transduction pathway activation studies; 4) reduction in tissue depletion; and 5) decreased reagent utilization. Multiplex immunofluorescence is not widely used as an analytical technique for clinical samples due to its lower sensitivity than chromogenic methods and concerns about autofluorescence in formalin-fixed tissue. Through the use of multispectral imaging, the absorption pattern of overlapping chromogens can be spectrally unmixed into quantitative individual dye components. This process works accurately when chromogens are developed to a lower intensity than typical for most standard pathology laboratories. The chromogens must be developed to a point where they do not scatter light or block subsequent chromogenic development. The current study was performed to define the accuracy of spectral unmixing of chromogens and evaluate whether interference existed between chromogens. Stained samples were imaged using multispectral imaging, and quantitative grayscale images for each analyte were created by spectral unmixing. The grayscale images were analyzed using image analysis software to determine the overall expression levels of each analyte. Development with Vulcan Red followed by DAB produced little interference, while high density development with DAB decreased the subsequent Vulcan Red staining intensity. Vulcan Red, DAB and hematoxylin unmixed properly and quantitatively. In summary, when multiplexed chromogenic IHC assays were developed using methods optimized for multispectral imaging, the technique demonstrated excellent performance characteristics without suffering from drawbacks commonly associated with fluorescence-based techniques.

Introduction and Methods

Multiplex immunohistochemistry (IHC) analysis was performed with Ki-67 and an antibody to Antigen #2 (target name withheld). Multiplex immunohistochemistry was performed using a cocktail of primary antibodies and polymer-based detection chemistry. Staining for Ki-67 was developed first with alkaline phosphatase and Vulcan Red chromogen. Staining for Antigen #2 was developed second using HRP and DAB. Samples were counter-stained with hematoxylin. If

Spectral Imaging


1. Pathologist reads this image (DAB + Vulcan Red, Ki-67 + DAB, Hematoxylin counterstain).

2. Using spectral unmixing, image stack is unrolled into grayscale images for each stain.


3. Composite images are re-colored for ease of interpretation (Ki-67=Red, Antigen #2=DAB, Hematoxylin (Htx) counterstain).

deposition of Vulcan Red masked the development of the second stain, there would be a reduction in the intensity of DAB where it co-localized with Vulcan Red, creating a "donut-like" appearance. If DAB caused inappropriate unmixing or masking of the second stain, the Ki-67 staining intensity would be reduced in cells that were double-stained with DAB for Antigen #2 compared to cells that were only Ki-67 positive.

Results



Ki-67 immunofluorescence using Alexa 488 labels basal layer of squamous epithelium (Sq) and germinal center B cells (G) in human tonsil. Note the autofluorescence (A) in formalin-fixed paraffin-embedded tissue.




Strong green autofluorescence due to collagen and light green general autofluorescence in formalin-fixed paraffin-embedded tissue can interfere with analysis.

Comparison of Immunofluorescence to Chromogenic Multiplex IHC

	Fluorescent	Chromogenic
Stability	Dyes Photobleach and degrade over time Quantum dots require "charging up" on fluorescent lightbox to reach "steady" intensity	Stable over time
Convenience	Requires imaging in the dark	Brightfield microscopy
Sensitivity	Limited or no signal amplification Large quantum dots = steric hindrance (esp. nuclear stains)	Enzymatic detection provides signal amplification
Endogenous Stain/Fluorescence	Some autofluorescence can be subtracted out using software, but it is not entirely spectrally uniform (e.g. collagen, elastin, lipofuscin)	Melanin in skin and other tissues can interfere with DAB-based detection
Number of Analytes	Theoretically, can perform 6-plex or greater	Have unmixed 4 chromogens successfully
Linearity	Linearity has been demonstrated	Linearity has been demonstrated


Unmixed Multiplex Chromogenic Images



Ki-67=Vulcan Red, Antigen#2 (An#2)=DAB, Hematoxylin (Htx) counterstain.

Multispectral cube taken and unmixed into components. Re-colored Ki-67=Red, An#2=Green, Htx=blue

Lack of Interference by First Round Stain




Note the uniformity of staining for Antigen#2 (DAB), despite the underlying Ki-67 (Red) staining.

If the first round chromogenic deposition interfered with the development of the second round chromogen, the green would have lighter intensity in the center where it overlapped the red.

There is no evidence that the deposition of Vulcan Red reduces the staining for the second analyte.

Lack of Interference by Second Round Stain



Note the uniformity of staining for Antigen#2 (DAB), despite the underlying Ki-67 (Red) staining.

The staining for Ki-67 (Red) was evaluated in regions containing a high fraction of DAB co-localization vs. regions with Ki-67 alone (demonstrated above).

Although variability exists between regions, there was no evidence that DAB masked the underlying red stain.

Region (ROI)	0631-0728	0281-0938	0675-0479	0122-0417
Description	NO DAB	NO DAB	DAB	DAB
Mean	14.21	20.05	17.05	22.03
StDev	19.08	22.96	19.32	25.61
Min	0	0	0	0
Max	111	115	103	113
Mode	0	0	0	0
Total	31396	29483	40877	45098
OD 20-39	13.55%	17.66%	18.60%	15.11%
OD 40-59	7.33%	9.21%	8.82%	10.48%
OD 60+	10.20%	15.75%	12.12%	20.05%

Conclusions

- Development with Vulcan Red followed by DAB produced little interference.
- Vulcan Red, DAB and hematoxylin unmixed properly and quantitatively.
- Multiplexed chromogenic IHC assays were developed using methods optimized for multispectral imaging. The technique demonstrated excellent performance characteristics without suffering from drawbacks commonly associated with fluorescence-based techniques.